

<b>2025 Blood and Urine Drug Detection Levels for Randox Screening</b>			
<b>Compound</b>	<b>Blood (ng/mL)</b>	<b>Urine Tox (ng/mL)</b>	<b>Urine OWI* (ng/mL)</b>
Cocaine Metab.	50	50	<b>188</b>
Opiates	10	50	<b>2817</b>
Amphetamine	20	500	<b>582</b>
Methamphetamine	20	500	<b>685</b>
Amphetamines Class**			<b>685</b>
Cannabinoids	10	50	<b>61</b>
Benzodiazepine 1	20	50	50
Benzodiazepine 2	10	50	50
Benzodiazepine 3	10	50	50
PCP	10	25	<b>38</b>
Methadone	50	150	150
Barbiturates	200	200	200
Tricyclics	50	100	100
Buprenorphine	10	5	5
Oxycodone 1	10	50	<b>233</b>
Oxycodone 2	10	50	<b>233</b>
Gen Opioids	10	100	<b>385</b>
Dextromethorphan	20	25	25
Meprobamate	500	250	250
Zolpidem	10	10	10
Tramadol	10	10	10
Fentanyl	2	2	2

\*Urine OWI levels include margin of error

\*\*Only utilized in Urine OWI cases